





Objectives:

This experiment aims to:

- I. determine the Michaelis-Menten constant (K_m) and the maximum rate of reaction done with tyrosinase (V_{max}) via spectrophotometry
- II. determine the type of inhibition done (competitive, non-competitive, or mixed) by p-hydroxybenzoic acid and a selected household inhibitor on tyrosinase


Reagents:

<i>Reagent Name and Structure</i>	<i>Physical and Chemical Properties</i>	<i>Safety Precautions</i>
4.80 mM catechol stock solution 	Brown powder MP: 100°C Density: 1.34 g/cm ³ at 15°C MW = 110.11 g/mol Absorption maximum at 280 nm; does not absorb significantly at 410 nm	Causes serious skin and eye irritation if in contact Wear protective gears (i.e. protective gloves, lab gown) Rinse immediately with water if in contact Disposal: G704
0.10 M pH 7 Phosphate (H ₃ PO ₄ and NaH ₂ PO ₄) 	- odorless liquid - density: 1.91 g/cm ³	Rinse contaminated skin and eyes with running water Disposal: B204 (if excess is <pH7)
Tyrosinase 	Solid state MP: 89-94°C catalyzes the oxidation of catechol to 1,2-benzoquinone, of which absorbs 410 nm strongly	Rinse eye and skin with running water if contaminated Avoid inhaling the compound Disposal: M503
0.30 mM p-hydroxybenzoic acid 	Beige powder solid Odorless MP: 214°C MW = 138.12 g/mol	Hazardous combustion compounds Rinse cautiously with water if ever the skin and eyes are contaminated Disposal: B208


Deionized water $\begin{array}{c} \text{H} & & \text{H} \\ & \diagdown & / \\ & \text{O} \end{array}$	- colorless, odorless liquid - BP: 100°C - MW: 18.02 g/mol	Non-irritant, generally safe to handle Disposal: sink

Procedures:

Part A. DETERMINATION OF K_m and V_{max}



Part B. ENZYME INHIBITION STUDIES



Calculations:

I. Calculations for the Preparation of Catechol Stock Solutions

To prepare 2.40, 1.20, 0.60, and 0.30 mM catechol solutions from 4.80 mM stock catechol, given that each solution is limited to 15 mL:

- For 2.40 mM solution (A): 7.5 mL stock + 7.5 mL DI H₂O
- For 1.20 mM solution (B): 7.5 mL soln. A + 7.5 mL DI H₂O
- For 0.60 mM solution (C): 7.5 mL soln. B + 7.5 mL DI H₂O
- For 0.30 mM solution (D): 7.5 mL soln. C + 7.5 mL DI H₂O

These values are obtained from the dilution equation:

$$M_1V_1 = M_2V_2$$

Example for 2.40 mM solution (A):

$$V_1 = M_2V_2/M_1 = (2.40 \text{ mM})(15 \text{ mL}) / (4.80 \text{ mM}) = 7.5 \text{ mL Stock}$$

II. Calculation of Molar Concentration for Each Batch of Solutions

Assuming solution volumes are additive, the volumes of each standard solution in the 96-well plate are:

- For tyrosinase-only batch: $V = 150 \mu\text{L} + 50 \mu\text{L} = \underline{200 \mu\text{L}}$
- For tyrosinase-and-inhibitor batch: $V = 150 \mu\text{L} + 50 \mu\text{L} + 50 \mu\text{L} = \underline{250 \mu\text{L}}$

Because 150 μL of the catechol solution is added per well, and that the solution is diluted to 200 or 250 μL , the concentration of catechol is decreased, also in accordance to the dilution equation:

$$M_1V_1 = M_2V_2$$

Example for 2.40 mM solution (A) in Tyrosinase-and-Inhibitor Batch:

$$M_2 = M_1V_1/V_2 = (2.40 \text{ mM})(150 \mu\text{L}) / (250 \mu\text{L}) = 1.44 \text{ mM Catechol}$$

III. Notes and Formulas for Data Analysis

- The Beer-Lambert law can be used to quantify the molar concentration c of catechol before the reaction and 1,2-benzoquinone after the reaction. Since the absorbance A at 410 nm is measured, and 1,2-benzoquinone, but not catechol, absorbs light at 410 nm (large molar absorptivity ϵ), the reacted tyrosinase-only solutions obtained after 7 minutes can be set as the external standard solutions from which an external calibration curve may be obtained. It is assumed here that all the catechol is transformed to 1,2-benzoquinone and that the original concentration of catechol is equal to the final concentration of 1,2-benzoquinone. Note that b is path length through each well and A_0 is systematic error.

$$A = A_0 + \epsilon bc$$

- Using the above equation, the molar concentrations of 1,2-benzoquinone for each well at each time value may be calculated. Once the concentration is plotted as y against time t for a particular solution (tyrosinase with or without inhibitor), the points may be fitted with an exponential model (in the form $y = A - Ae^{-kt}$) and the reaction velocity v may be calculated as the magnitude of the slope, $|dy/dt|$, as time x approaches zero.

$$v = \left| \lim_{t \rightarrow 0} \left[\frac{dy}{dt} \right] \right|$$

- The initial molar concentration of substrate, $[S]_0$ is equal to the initial concentration of catechol in each well, which is calculated from the concentrations of standards used based on the volume of the solution (see calculations II and table 3).
- Values of $1/v$ and $1/[S]_0$ are plotted against each other as y and x axes, respectively. This is called a Lineweaver-Burk plot. Then, the maximum velocity v_{max} and Michaelis-Menten constant K_M may be calculated according to the Lineweaver-Burk equation. The first equation is used for the tyrosinase-only samples, whereas the second equation is used for samples with tyrosinase and inhibitor (Atkins & de Paula, 2006)

$$\frac{1}{v} = \frac{1}{v_{\max}} + \left(\frac{K_M}{v_{\max}} \right) \frac{1}{[S]_0}$$

$$\frac{1}{v} = \frac{\alpha'}{v_{\max}} + \left(\frac{\alpha K_M}{v_{\max}} \right) \frac{1}{[S]_0}$$

- Then, the $1/v$ vs $1/[S]_0$ plots of tyrosinase-only (TO) and tyrosinase-and-inhibitor (TI) samples may be compared. The type of inhibition may be identified as follows (Atkins & de Paula, 2006):
 - If the slope of plot TI is significantly greater than the slope of plot TO, the inhibitor is a **competitive inhibitor** for tyrosinase.
 - If the intercept of plot TI is significantly greater than the intercept of plot TO, the inhibitor is a **non-competitive inhibitor** for tyrosinase.
 - Finally, if both the slope and intercept of plot TI are significantly greater than the slope and intercept of plot TO, the inhibitor is a **mixed inhibitor** for tyrosinase.

Gantt Chart:

Each portion must be done after the previous step is finished; Each portion has its own allocated time

[illegible]

Insert Solutions + Inhibitor to Well Plates																
Travel Time																
Addition of Tyrosinase																
Spectrophotometry																

Data and Observations:

- Setup of the 96-Well Plate (indicates where each solution is located in)

Table 1. *Measurements Taken During Preparation of Catechol Solutions*

Preparation of Catechol Stock Solution	
Parameter	Measurement
Solution pH	

Table 2. *Setup of The 96-Well Plate*

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

Table 3. *Legend of Catechol-containing External Standard Solutions to be Prepared*

Aspect	Label	Solution A	Solution B	Solution C	Solution D
--------	-------	------------	------------	------------	------------

Catechol	150 μ L mM Catechol Added	2.40	1.20	0.60	0.30
	Initial mM Catechol for Tyrosinase-Only Batch	1.80	0.90	0.45	0.225
	Initial mM Catechol for Tyrosinase-and-Inhibitor Batch	1.44	0.72	0.36	0.18
Total Volume of Solution	Total Volume (μ L) for Tyrosinase-only Batch	200	200	200	200
	Total Volume (μ L) for Tyrosinase-and-Inhibitor Batch	250	250	250	250

Table 4. Absorbance of Catechol-1,2-Benzoquinone Reaction Mixtures at 410 nm for Tyrosinase-Only Batch

	Absorbance at 410 nm (Tyrosinase Only)			
Elapsed Time (s)	Solution A	Solution B	Solution C	Solution D
0				
15				
30				
45				
60				
:	:	:	:	:
420				

Table 4. Absorbance of Catechol-1,2-Benzoquinone Reaction Mixtures at 410 nm in the Presence of Tyrosinase and 0.30 M *p*-Hydroxybenzoic Acid

	Absorbance at 410 nm (Tyrosinase and <i>p</i> -HBA)			
Elapsed Time (s)	Solution A*	Solution B*	Solution C*	Solution D*
0				
16				
32				

48				
64				
⋮	⋮	⋮	⋮	⋮
432				

Table 5. Absorbance of Catechol-1,2-Benzoquinone Reaction Mixtures at 410 nm in the Presence of Tyrosinase and a Household Inhibitor

	Absorbance at 410 nm (Tyrosinase and _____)			
<i>Elapsed Time (s)</i>	<i>Solution A'</i>	<i>Solution B'</i>	<i>Solution C'</i>	<i>Solution D'</i>
0				
16				
32				
48				
64				
⋮	⋮	⋮	⋮	⋮
432				

Observations:

References:

- Albarrán, G., Boggess, W., Rassolov, V., & Schuler, R. H. (2010, June 25). Absorption spectrum, mass spectrometric properties, and electronic structure of 1,2-benzoquinone. *Journal of Physical Chemistry*. <https://pubs.acs.org/doi/10.1021/jp101723s>
- Atkins, P., and de Paula, J. (2006). *Atkins' Physical Chemistry*. W. H. Freeman and Company.
- Deionized Water. (n.d.). Fisher Scientific. Retrieved September 17, 2024, from <https://fscimage.fishersci.com/msds/89955.htm>
- Safety Data Sheet - 4-Hydrobenzoic Acid. (2010, November 8). Thermo Fisher Scientific Safety Data Sheet. Retrieved October 22, 2024, from <https://www.fishersci.com/store/msds?partNumber=AC120991000&productDescription=4-HYDROXYBENZOIC+ACID%2C+9+100GR&vendorId=VN00032119&countryCode=US&language=en>
- Scalzone, A., Bonifacio, M. A., Cometa, S., Cucinotta, F., de Giglio, E., Ferreira, A. M., & Gentile, P. (2020, June). pH-Triggered Adhesiveness and Cohesiveness of Chondroitin Sulfate-Catechol Biopolymer for Biomedical Applications. *Frontiers in Bioengineering and Biotechnology*. 10.3389/fbioe.2020.00712
- Sigma Aldrich. (n.d.). *Pyrocatechol - Safety Data Sheet*. Sigma Aldrich. Retrieved October 15, 2024, from <https://www.sigmaaldrich.com/PH/en/sds/SIAL/C9510?userType=undefined>
- Phosphate Buffered Saline. (2022). Sigma-Aldrich. Retrieved September 17, 2024, from <https://www.sigmaaldrich.com/PH/en/sds/SIGMA/P5119?userType=undefined>
- Tyrosinase from Mushroom. (2023). Sigma-Aldrich. Retrieved September 17, 2024, from <https://www.sigmaaldrich.com/PH/en/sds/sigma/t3824?userType=undefined>